





Housing Horn Angles



0° 15° 45° 90°



135° 180° 225° 270°

B-155 B-165 B-185
VIEW FROM ANODE END OF THE TUBE

Euro 270° Euro 90°

Toshiba Housings
VIEW FROM ANODE END OF THE TUBE

X GX FX

Coverage Chart

Focal to Film Distance	TARGET ANGLE			
	7°	10°	12°	16°
(in)				
100	24.6	35.3	42.5	57.4
80	19.7	28.2	34	45.9
72	17.7	25.4	30.6	41.3
60	14.7	21.2	25.5	34.4
48	11.8	16.9	20.4	27.5
44	10.8	15.5	18.7	25.2
42	10.3	14.8	17.9	24.1
40	9.8	14.1	17	22.9
38	9.3	13.4	16.2	21.8
34	8.4	12	14.5	19.5
30	7.4	10.6	12.8	17.2
28	6.8	9.7	11.7	15.8
26	6.3	9	10.9	14.7
24	5.8	8.3	10	13.6